

Optimization of IPRP-UPLC Method to Prevent PS to PO Conversion on Column

Jiaqi 'Jack' Zheng¹, Derek Gauntlett¹, Asaf Cohen² and Gili Hart²

¹LGC Biosearch Technologies, 2199 South McDowell Blvd, Petaluma CA 94954, USA

²Biohouse Labs, Minrav Bldg. 1st Floor, Haddasah Ein Karem, Jerusalem, Israel 9112101


Jiaqi.zheng@lgcgroup.com

Abstract

Several chemical modifications have been developed in antisense therapeutic drugs to date. The phosphorothioate (PS) bond is one of the original and most widely-used backbone variants. PS bond oligo-modification alters the phosphodiester (PO) bond by replacing one of the non-bridging

oxygens with a sulfur atom. This modification improves the metabolic stability and cellular uptake of oligonucleotide without affecting their affinity for target. During development of the analytical IPRP method for compound X, which contains PS backbone, we observed a phenomenon that PS bond was being

oxidized to a PO bond on-column. Here, we present an approach to improve an IPRP method that prevents on column PS to PO oxidation by utilizing synergistic effects of temperature, addition of EDTA and DTT into the mobile phase.



Materials & Methods

Instruments:

- Waters ACQUITY UPLC H-Class
- Xevo G2-S ToF



Column:

- Waters ACQUITY Premier Oligonucleotide C18 Column, 130Å, 1.7 µm, 2.1 x 100 mm, Part # 186009485

Reagents:

- Water, LC/MS Grade, Fisher, W6-4
- HFIP, ≥ 99%, TCI, H0424
- DIPEA, 99.5%, Sigma, 496219
- Acetonitrile, LC/MS Grade, Fisher, A955-1
- DTT, 98%, Sigma, D9779
- EDTA, ≥99.5%, Fisher, #BP118-500

Compound X:

- Oligonucleotide with PS bond modification

Software:







- Empower 3
- MassLynx

Original Method Condition:

- Column Temp: 80 °C
- Mobile Phases:
 - A: 1% HFIP, 0.1% DIPEA in water
 - B: 40% Acetonitrile in water

Time	%A	%B	Curve
Initial	98.0	2.0	*
2.0	98.0	2.0	Linear
10.0	71.0	29.0	Linear
32.0	61.0	39.0	Linear
33.0	50.0	50.0	Linear
34.0	98.0	2.0	Linear
41.0	98.0	2.0	Linear

Results



Final Method

- Column Temp: 60 °C
- Mobile Phases A: 1% HFIP, 0.1% DIPEA, 10 µM of EDTA, 100 µM of DTT in water.
- Mobile Phases B: 40% Acetonitrile, 10 µM of EDTA, 100 µM of DTT in water.
- LG Gradient:

Time	%A	%B	Curve
Initial	98.0	2.0	*
2.0	98.0	2.0	Linear
10.0	69.0	31.0	Linear
32.0	59.0	41.0	Linear
33.0	40.0	60.0	Linear
34.0	40.0	60.0	Linear
35.0	98.0	2.0	Linear
42.0	98.0	2.0	Linear




Fig. 8: LC chromatogram of standard solution injection using the final method. The FLP purity is ~92% and the %Area for PS to PO impurity is ~2%.




Fig. 9: Total Ion Count TOF mass spectrum of standard solution injection using the final method. Addition of DTT and EDTA in the mobile phases did not affect the MS responses. The MS signal for the FLP peak is as high as 5.53 X 10⁶.

Conclusion

- The PS bond on compound X appears to favor being converted to PO using the originally developed method.
- The PS to PO oxidation is probably attributed by the high temperature and cation metals.
- Lower temperature, metal chelating agent (EDTA) and reducing agent (DTT) all contribute to prevent the PS to PO conversion.
- Addition of DTT and EDTA in the mobile phases do not affect the UV and MS responses.